

**Bovine Serum Albumin**  
**Standard Grade Powder ph. 7.0**  
**Manufactured in New Zealand**

Proliant NZ Ltd Bovine Serum Albumin (BSA) is a high purity, lyophilized "Fraction V" powder derived from the plasma of MPI-inspected healthy New Zealand origin animals. The albumin is isolated from other plasma protein products and lipids by a unique proprietary, non-solvent procedure. The separation process includes a modified Cohn / heat-shock procedure designed to inactivate proteases and other potentially interfering enzymes. Further processing such as extensive membrane dialysis and filtration, minimizes analytes that can cause background interference in sensitive in vitro diagnostic (IVD) assays as well as inhibition in sensitive cell and microbial culture systems. Typical purity levels are in excess of 99%. Proliant's New Zealand Source Standard Grade bovine albumin product is an economical choice that is ideal for most diagnostic and research applications.

**APPLICATIONS**

- Protease-sensitive Immunoassays such as RIA, EIA, Fluorescent and Chemiluminescent
- Protein Standard, Diluent
- Protein, Conjugate or Enzyme Stabilizer
- Hybridization
- Selected Cell Culture applications

**FEATURES / BENEFITS**

- High Purity, Low Interference, Low Background
- Virtually No Detectable Protease Activity, Ensures Assay Integrity
- "Closed Loop" Manufacturing\*
  - Minimizes Contamination and Maximizes Reproducibility
- Superior Solubility/Filterability for Ease of Use
- Product Development and Technical Services Designed for maximum responsiveness and Flexibility
- Manufactured from MPI inspected Bovine Plasma
- Manufactured According to FDA cGMPs
- Validated for Prion (TSE) Clearance

**SCIENTIFIC ANALYSIS**

	<b><u>Specification</u></b>	<b><u>Typical Values</u></b>
Protein (dry basis)	≥ 96.0%	99%
Purity (albumin)	≥ 98%	100%
Solubility (4% water)	Clear/slightly hazy	Clear/slightly hazy
Moisture	≤ 5.0%	< 2%
pH (10% water)	6.5-7.5	6.8-7.2
Ash	< 2%	<2%
Protease	≤ 0.005 Units/mg	≤ 0.005 Units/mg
Endotoxin	≤ 3 EU/mg	0.1 EU/mg
IgG	Not Detected	Not Detected

**PHYSICAL CHARACTERISTICS**

Proliant NZ Bovine Serum Albumin is a highly soluble yellow to green with tan to green cast lyophilized powder. A solution of up to 30% in de-ionized water is clear to slightly hazy and virtually particulate-free.

**PACKAGING**

Standard package size is 10kg.  
Smaller sample sizes for evaluation are available.

**STORAGE RECOMMENDATIONS**

Store sealed in a cool, dry environment for 3 years, and up to 5 years with re-qualification.

**ORDERING INFORMATION**

Please order by product name "Bovine Serum Albumin, New Zealand Standard Grade pH 7.0 Lyophilized Powder" and refer to SKU # 69100. Proliant New Zealand Limited 27-33 Kawakawa Road, Feilding, 4775, NEW ZEALAND PO Box 366, Feilding, 4740, NEW ZEALAND. Phone: 64 6 280 4150. [www.proliantinc.com](http://www.proliantinc.com)

\*Proliant's unique "Closed Loop" Manufacturing means that we collect only in exclusive abattoirs with our own approved proprietary collection system. The material undergoes ante-and post-mortem inspection to ensure it is the highest edible grade before it is shipped to our manufacturing facility. No other manufacturer can claim this important closed link, as they must buy their raw material from outside sources and the spot market.

What's more, from the time our product enters our tanks to the moment it is loaded into the dryers, the product is not exposed to the external environment. Once the raw material reaches our cGMP plant, the material is loaded directly into our process and manufactured through a series of closed tanks, lines, filters and separation equipment. In our plant, there are no exposed filter presses or open tanks, nor is any non-food grade or noxious solvent material introduced into the process. There is also no possibility of contamination with other animal-derived products or bovine tissues in our completely dedicated facility.